




# Solutions for the Production of Solids

One Partner. Maximum Efficiency.



**FETTE**  
COMPACTING

# together

## from lab to production

### Our Strategic Approach

- 4 together – from lab to production**  
Faster time-to-market with a holistic process partnership
- 6 Quality by Design and Qualified Experts-Database**  
The foundation for efficient production of solids

### Research and Development

- 8 R&D Solutions: The Reliable Path to Market Readiness**  
Data-driven. Scalable. Safe.
- 9 Powder Compaction Analyzing Unit F Lab 10**  
Understanding powders before production begins
- 10 Emulators**  
Process development close to production
- 12 Galenics-Tablet Press 102i**  
Production-related galenics for direct scale-up
- 14 R&D Services**  
Our formula for your efficient research

### Better together – from lab to production

For more than 75 years, Fette Compacting has been redefining precision in industrial tablet production. As the world market leader in integrated solutions for OSD (Oral Solid Dosage) manufacturing, we deliver decisive competitive advantages. Our portfolio ranges from high-performance tablet presses to specialized Containmentment solutions and process equipment.

But technology is only one part of our DNA. Guided by our principle “Together – from lab to production”, we support you throughout the entire life cycle of your products: from initial formulation in the lab through technology transfer to commercial production.

That way, you find the right solution for every product and get it to market faster. With our global network of five Competence Centers in Germany, China, the USA, India, and Brazil, along with 13 subsidiaries and representatives in 50 countries, we are always right where you need us.

As your holistic process partner, we provide complete solutions. Whether pharma, nutrition, or chemical: We combine tailored consulting with deep technological know-how. The result: customized, efficient solutions for every phase of your product life cycle – and a partnership that delivers real value.

### Production

- 16 Services for Production**  
Close to your machine – just like you
- 18 FE CPS**  
The new standard in Continuous Direct Compression
- 20 ePAT**  
Integrated process analytics for Continuous Manufacturing and batch production
- 22 FE Series**  
Flexibility at the highest level
- 23 FE35
- 24 FE55
- 25 FE75
- 26 i Series**  
Versatile options, no compromises
- 27 F10i
- 28 F20i
- 29 F30i
- 30 p Series**  
Simply efficient production
- 31 F10p
- 32 F20p
- 33 F30p

- 34 Containment Guard**  
Processing active ingredients safely
- 35 FE55 with Containment Package
- 36 FE75 with Containment Package
- 37 F10i Containment Guard
- 38 F20i Containment Guard
- 39 F30i Containment Guard
- 40 Wash in Place**  
Clean safely, produce faster
- 41 1090i WiP
- 42 2090i WiP
- 43 3090i WiP
- 44 Tableting Tools**  
Always in top form
- 46 Process Equipment**  
Seamless processes in production and active ingredient analysis
- 48 Handling Systems**  
Fast changeovers – keep things moving
- 49 EasyCare**  
Safe and efficient handling of tableting tools
- 50 SmartInterface**  
Real-time data, anytime and anywhere



### General

- 51 Global Training Program**  
Knowledge that works. From lab to production.
- 52 Next Level Nutrition**  
Tableting solutions for food and supplements

# Together – from lab to production

Faster time-to-market with a holistic process partnership

Research and Development



Extended range of services along the process chain:



## Knowledge Base QED

- › Concentrated expertise: Access to more than 75 years of tableting experience
- › Confident decisions: The digital foundation for your well-informed process planning


## Formulation Development

- › Formulation consulting: Expert support backed by pharmaceutical expertise
- › Understand your materials: Precise characterization of your active ingredients and excipients

## Process Development

- › Test without risk: Emulator-based trials
- › Stay in control: Development of tailored control strategies
- › Leverage your data: Seamless integration of Process Analytical Technology (PAT)

Production



## Technology Transfer

- › Scale with ease: Reliable transition from development to production
- › Navigate regulatory requirements: Support in overcoming regulatory hurdles
- › Validate your processes: Efficient transfer and safeguarding of your processes

## Production

- › Produce flexibly: Tablet presses for batch operations and Continuous Manufacturing
- › Protect people and the environment: Purpose-built Containment solutions
- › All from a single source: Process equipment and tableting tools

## Process Optimization

- › Boost performance: Ongoing production support
- › Minimize risks: Identify and address weak points early
- › Improve processes: Development of customized optimization strategies



Brochure: Together – from lab to production

# Quality by Design and Qualified Experts-Database

The foundation for efficient production of solids


## Quality is not a coincidence — it is built into the design

Excellent tablets are the result of systematic planning, which is why we embrace Quality by Design (QbD). By understanding all key production elements and material properties, we ensure that quality is not just checked at the end of production, but is embedded in the product from the very beginning..

For robust processes, we examine the interactions between Critical Material Attributes (CMAs), Critical Process Parameters (CPPs), and Critical Quality Attributes (CQAs). By understanding how material and machine interact and how they influence quality, we work with you to define the proven process space where your production runs consistently, reproducibly, and at flawless quality.



## Optimized formulation and process development with QED



QTPP	CQAs	CMAs	CPPs	CPPs ✓
Customer provides materials and requirements		Material characterization and analysis	Prediction of process parameters	Verification and optimization of predicted process parameters

QTPP = Quality Target Product Profile  
 CQAs = Critical Quality Attributes  
 CMAs = Critical Material Attributes  
 CPPs = Critical Process Parameters

**The optimal solution with minimal investment of time and material**

### A knowledge edge at every stage of development

Quality by Design delivers the plan; QED delivers the answers. With the Qualified Experts-Database (QED), Fette Compacting has built an extensive data platform that draws on more than 75 years of tableting experience. It links formulation and process expertise and serves as a solid decision-making basis throughout development and production. AI-driven analyses enable early predictions for optimal process development. Our experts verify these predictions and advise you based on the QED with maximum efficiency. This saves you valuable time on the path to market readiness.

### The strengths of the QED:

- + Links multivariate material data with real production and quality data
- + Covers an exceptional range of machines, format parts, and process variants
- + Uses AI to reliably analyze and predict complex interactions
- + Significantly reduces development effort, trial runs, and material consumption
- + Delivers concrete, practice-oriented recommendations for robust and efficient processes

# R&D Solutions: The Reliable Path to Market Readiness

Data-driven. Scalable. Safe.

The path from the first formulation idea to industrial production is complex and often fraught with risk. Critical questions arise: Will the powder behave on the high-performance machine the same way as in the lab? Are the processes stable and economically scalable?

With Fette Compacting's R&D Solutions, you significantly accelerate this process. We offer a seamless system of laboratory equipment, process expertise, and services to eliminate uncertainties early on.

Instead of catching deviations during ongoing production, we identify them in the development phase. Whether through precise lab analyses, emulator trials, or test runs on the galenics-tablet press: We provide the physical and data-driven foundation for informed decisions. The result: robust processes and a significantly shorter time-to-market.

# Powder Compaction Analyzing Unit F Lab 10

Understanding powders before production begins

To avoid costly errors in tableting, full transparency into material properties is essential. The F Lab 10 characterizes the compaction behavior of powders with minimal material usage.

The system comprehensively analyzes active ingredients and excipients and delivers precise data for the entire product life cycle. It enables early identification of critical material attributes and the development of robust formulations before production starts on large-scale machines.



## **F**<sub>Lab 10</sub>

<b>Maximum load</b>	1,000 kg
<b>Punch size range (diameter)</b>	3–15 mm
<b>Compaction speed</b>	0.01–3 mm/s
<b>Data capture rate</b>	200 Hz
<b>Load cell travel</b>	40 mm
<b>Load cell resolution</b>	1:5,000
<b>Calibration</b>	Dead weights or proving ring
<b>Power requirement</b>	90–240 VAC 3.15 A
<b>Dimensions</b>	320 × 285 × 388 mm
<b>Weight</b>	24 kg






Brochure: R&D Solutions

# Emulators

## Process development close to production

The scale-up from lab to production is a critical phase. Simulators reach their limits here because they can only replicate complex physical effects to a limited extent. Our emulators, by contrast, faithfully reproduce the physics and mechanics of production-scale machines – without tying up valuable capacity on your high-performance equipment.

**The advantages:** You reliably validate critical parameters and analytical methods, save expensive active ingredients, minimize scale-up risks, and transfer the data directly to the production line.



### FE CPS Process Emulator

#### Precision in Continuous Dosing

The FE CPS Process Emulator physically replicates the dosing process of the production machine (FE CPS) at a 1:1 scale. It uses identical Loss-in-weight feeders and refill systems. This allows dosing parameters, screw types, and refill strategies to be systematically tested.

- + Direct scale-up: Recipe transfer via “plug-and-play” without conversion factors
- + Independence: Test refill strategies and flow behavior without tying up ongoing production
- + Material efficiency: Process optimization with minimal raw material usage

<b>Electrical supply parameters</b>	Operating voltage 400–480 V, Frequency 50–60 Hz
<b>Power supply</b>	3 phase + PE
<b>Dimensions L x W x H</b>	1,642 x 822 x 2,047 mm
<b>Weight</b>	650 kg
<b>Supply air volume flow</b>	50 m <sup>3</sup> /h
<b>Compressed air supply</b>	6 bar, 300 l/min
<b>Exhaust air volume flow</b>	50 m <sup>3</sup> /h
<b>Connection to exhaust air system</b>	33.7 mm
<b>Network connection</b>	RJ45
<b>Display</b>	19" touch display
<b>Number of dosing stations</b>	1
<b>Throughput capacity</b>	0.7–340.1 kg/h*

\* depends on configuration

### ePAT Emulator

#### The Foundation for Reliable In-line Measurements

This emulator establishes the foundation for controlling tablet uniformity within the embedded Process Analytical Technology (ePAT). The system precisely emulates the dynamics of the tablet press: Tablets pass the sensor at the same speed and trajectory as in the actual process.

- + Process reliability: Validation of analytical models before production starts
- + Precise analytics: Optimization of spectrometer settings under real dynamic conditions
- + Resource efficiency: Representative tests and calibrations with a minimal number of tablets

<b>Electrical supply parameters</b>	Operating voltage 230 V, Frequency 50 Hz
<b>Dimensions L x W x H</b>	1,656 x 592 x 1,027 mm
<b>Weight</b>	298–337 kg*
<b>Main air inlet/outlet</b>	98 m <sup>3</sup> /h
<b>Air inlet/outlet</b>	19 m <sup>3</sup> /h
<b>Network connection</b>	RJ45
<b>WLAN</b>	Current and older WiFi standards
<b>Display</b>	19" touch display

\* depends on configuration

### ePAT BU Emulator

#### Understanding Powder Dynamics in the Feeding System

The system physically replicates the feeding device and powder supply. An integrated BU sensor measures powder homogeneity directly in the feeding device. This enables powder behavior analysis without starting a compression cycle.

- + Deep process understanding: Development of spectroscopic methods under realistic filling conditions
- + Material efficiency: Resource-conserving reuse of uncompressed powder
- + Clear focus: Dedicated analysis of the filling process without tying up tablet press capacity


<b>Electrical supply parameters</b>	Operating voltage 400–480 V, Frequency 50–60 Hz
<b>Power supply</b>	3 phase + PE
<b>Dimensions L x W x H</b>	1,235 x 675 x 2,010 mm
<b>Weight</b>	325 kg
<b>Supply air volume flow</b>	150–200 m <sup>3</sup> /h
<b>Exhaust air volume flow</b>	150–200 m <sup>3</sup> /h
<b>Connection to exhaust air system</b>	47.5 mm
<b>Network connection</b>	RJ45
<b>Display</b>	19" touch display

 ePAT sensor position

# Galenics-Tablet Press 102i

## Production-related galenics for direct scale-up

The 102i tablet press is the ideal solution for lab operation and galenical development. All parameters determined in the lab can be directly transferred to production-scale machines.



### Tablet

Number of layers 1-, 2-, 3-layer  
 max. tablet output 230,400 units/h  
 max. tablet output Pmax® 324,000 units/h  
 max. tablet diameter 25 mm



Die (D) / Segments (S)		D	D	D	D	D	D	D	D	S	S	S	S	
<b>Number of punch stations</b>		6	6	16 (8+8)	16 (8+8)	20		32	30	24	21	24	30	45
<b>Punch type</b>		FS19® EU19 TSM19 B	EU1* TSM1* D	FS19® EU19 TSM19 B	EU1* TSM1* D	EU1* EU1*-441 TSM1* D		FS19® EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1*-441 EU1* TSM1*	FS19® EU19 TSM19	FS12®	
<b>Tablet output units/h</b>	min.	9,000	9,000	24,000 (12,000)	24,000 (12,000)	30,000		48,000	45,000	36,000	31,500	36,000	45,000	67,500
	max.	43,200	36,000	96,000 (48,000)	96,000 (48,000)	120,000		230,400	216,000	172,800	126,000	144,000	216,000	324,000
<b>Max. compression force 1*</b>	kN	80	80	80	80	80		80	80	80	80	80	80	34
<b>Max. compression force 2*</b>	kN	80	80	80	80	80		80	80	80	80	80	80	34
<b>Max. tablet diameter</b>	mm	16	25	18	25	25		11	13	18	25	25	18	11
<b>Max. filling depth 1st layer**</b>	mm	18	18	18	18	18		18	18	18	20	20	20	20
<b>Pitch circle diameter</b>	mm	280	280	280	280	280		280	280	280	280	280	280	280
<b>Turret rotation speed min.</b>	min <sup>-1</sup>	25	25	25	25	25		25	25	25	25	25	25	25
max. (laboratory operation)	min <sup>-1</sup>	120 (150)	100 (100)	100 (100)	100 (100)	100 (100)		120 (150)	120 (150)	120 (150)	100 (150)	100 (150)	120 (150)	120 (150)
<b>Die diameter</b>	mm	30.16	38.1	30.16	38.1	38.1		22	24	30.16	-	-	-	-
<b>Die-/segment height</b>	mm	22.225	23.8	22.225	23.8	23.8		22.225	22.225	22.225	25	25	25	25
<b>Punch shaft diameter</b>	mm	19	25.35	19	25.35	25.35		19	19	19	25.35	25.35	19	12
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6		133.6	133.6	133.6	133.6	133.6	133.6	133.6
Upper/lower punch		(133.35)	(133.35)	(133.35)	(133.35)	(133.35)		(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)
<b>Upper punch insertion depth</b>	mm	1-4 (8***)	1-4 (8***)	1-4 (8***)	1-4 (8***)	1-4 (8***)		1-4 (8***)	1-4 (8***)	1-4 (8***)	1-4 (8***)	1-4 (8***)	1-4 (8***)	1-4 (8***)
<b>Dimensions L × W × H</b>	mm	920 × 1,136 × 1,875												
<b>Weight</b>		Tablet press 1,700–2,500 kg, operating terminal 100 kg												
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, 50/60 Hz, power consumption 8.4 kW												

Theoretical values and technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly.

\* limited by punch properties

\*\* special filling depth available on request

\*\*\* multi-layer operation

# R&D Services

## Our formula for your efficient research

With Fette Compacting's R&D Services, you expand your capabilities in pharmaceutical development. Whether you face a specific challenge or want to optimize your entire development process: Our experts support you with in-depth know-how and cutting-edge infrastructure.

We flexibly support you wherever you need it: from the first material analysis through the development and optimization of formulations to reliable transfer into production. You benefit from fast results, proven strategies, and a collaborative partnership that frees up your resources and keeps your projects on track.



## R&D Services

Further information about our services and training offerings for R&D can be found in our "R&D Solutions" brochure.



Brochure: R&D Solutions

### Formulation Development

#### Material Characterization

We analyze physical and chemical raw material properties as the basis for scalable manufacturing processes.

#### Powder Compaction Analysis

Precise analysis of compaction behavior to optimize formulations and processes.

#### HPLC

Precise analytics for research, method development, validation, and quality control – including reference measurements for Process Analytical Technology (PAT).

### Process Development

#### Feasibility Studies

Early evaluation of formulations and processes to minimize risk before investment decisions.

#### Process Design Space

Defining the safe parameter range for stable and quality-compliant production processes.

#### State of Control

Development of models for process control and complete batch traceability.

#### Control Strategy

Defining optimal control parameters to ensure long-term production stability.

### Technology Transfer

#### Process Transfer

Safe, validated transfer of your processes from the lab to the GMP production environment.

#### Clinical Batch

Comprehensive support for GMP-compliant manufacturing and release of clinical trial batches.

#### ePAT Development

Implementation of data-driven in-line process controls using near-infrared spectroscopy (NIR).

# Services for Production

## Close to your machine – just like you

Our Production Services provide holistic support for your tablet presses and production lines – from seamless commissioning through ongoing maintenance to data-driven process optimization. Complex products and stringent regulatory requirements place the highest demands on operations, so you need a partner who understands your machines down to every detail. With world-wide support, fast spare parts supply, and predictive maintenance, we ensure that your machines operate at peak efficiency, in full regulatory compliance and with maximum cost-effectiveness. We make your production future-proof.



### Production Services

Further information about our services and training offerings for production can be found in our “Production Services” brochure.



Brochure: Production Services

#### Production

##### Ramp-up Support

We support you directly after the qualification phase at the machine to minimize start-up risks and ensure a fast, regulation-compliant production start.

##### Technical Customer Support

Our global expert network ensures maximum machine availability through fast remote support, on-site assistance, and scheduled maintenance.

##### Spare Parts Management

With our digital spare parts catalog (eCAT) and worldwide logistics network, we supply you quickly and accurately with original parts (OEM).

#### Process Optimization

##### FMECA Service

We proactively identify and assess potential risks to minimize downtime through optimized spare parts inventory.

##### Process Support Production

By systematically gathering and evaluating your production data, we uncover untapped potential and lay the groundwork for efficiency gains.

##### Fleet Management

We take a holistic view of your machine fleet and develop tailored strategies to optimize utilization and maximize your return on investment (ROI).

##### General Refurbishment

We offer professionally refurbished, individually configured used machines as a cost-efficient and high-performance alternative to new purchases – and on request, we also modernize and overhaul your existing machines.

##### Machine Modernization

With targeted performance upgrades and modernizations, we bring your proven machines up to the latest standards in technology and safety.

# FE CPS

## The new standard in Continuous Direct Compression




Pharmaceutical solid dosage production is undergoing a paradigm shift: The transition from classic batch processing to Continuous Manufacturing delivers considerable efficiency and quality gains. Until now, complex systems with large footprints have been a barrier to adoption. This is precisely where Fette Compacting steps in, rethinking Continuous Direct Compression from the ground up.

With the FE CPS, we introduce a ground-breaking solution based on a rigorous Quality by Design approach. Instead of large mixing tanks and complex transfer systems, the FE CPS integrates all process steps – from precise dosing through homogeneous mixing to reliable powder transport – into a single, standardized machine concept. Material flows continuously in a single stream and is synchronized and monitored in real time.

### Your Advantages at a Glance:

- + Compact and modular: The system fits into standard pharmaceutical production rooms without costly modifications and can be flexibly deployed as a stand-alone solution or integrated into existing lines.
- + Simple and fast: The proven TRI.EASY operating concept from Fette Compacting and tool-free disassembly enable product changeovers and cleaning in less than one shift.
- + Maximum safety: Dust-tight units and a strict separation of process and technical areas guarantee effective Containment.
- + Fully controllable: The fully integrated embedded Process Analytical Technology (ePAT) monitors blend homogeneity and tablet uniformity continuously in real time.



### FE/CPS

<b>Throughput capacity</b>	5 to 250 kg/h with a single rotary tablet press 5 to 400 kg/h with a double rotary tablet press (formulation and tablet dependent)
<b>Number of dosing stations</b>	1 to 6
<b>Number of blender mixing zones</b>	2
<b>Feeder-blender funnel configuration</b>	Various configurations – formulation dependent
<b>Ingredient characteristics</b>	Powder density: 0.2 ~ 1.0 kg/l *
<b>Ingredient inlets</b>	DN150 tri-clamp connection
<b>Electrical supply parameters</b>	Operating voltage 400 – 480 V, frequency 50/60 Hz ±5%
<b>Power supply</b>	3 phase + PE
<b>Power</b>	9 kW
<b>Compressed air supply</b>	6 bar / 300 l/min peak volume flow
<b>Air extraction rate</b>	50 m³/h
<b>Air extraction negative pressure</b>	-10 hPa
<b>Air extraction unit connection</b>	70 mm
<b>Dimensions</b>	Machine height: 2,300 mm, Ingredient inlet: 2,355 mm Total height including blend conveyor arm: 3,000 mm
<b>Weight</b>	approx. 4,900 – 5,200 kg **

Subject to technical changes.  
\* Assumption of typical values; others on request  
\*\* Depending on the exact configuration



Brochure: FE CPS

# ePAT

## Integrated process analytics for Continuous Manufacturing and batch production

Whether you produce continuously with the FE CPS or use classic batch processes: Our embedded Process Analytical Technology (ePAT) enables seamless in-line monitoring of critical quality parameters during ongoing production.

Near-infrared (NIR) sensors monitor powder blend homogeneity (Blend Uniformity) and Tablet Uniformity, detecting quality deviations in real time. This prevents the late rejection of entire batches, minimizes waste, and represents a decisive step toward Real-Time Release.

### Blend Uniformity


The BU sensor verifies the uniform distribution of active ingredients and excipients, performing up to ten independent measurements per second. The active measuring head can be flexibly positioned – for example at the mixer outlet, at the tablet press inlet, or directly in the filling device (Fill-O-Matic).

### Tablet Uniformity

The TU sensor inspects 100 percent of all tablets directly at ejection for uniformity. At speeds of up to 120 tablets per second, the system detects deviations and rejects defective individual tablets.




TU sensor



BU sensor

ePAT  
embedded Process  
Analytical Technology



F  
E | CPS



Brochure: FE CPS

# FE Series

## Flexibility at the highest level

The tablet presses in the FE Series deliver high performance, flexible applications, and process-reliable production. They ensure quality and efficiency in tableting while combining innovation with user-friendliness. The FE55 and the FE75 can be combined with the FE CPS for continuous production and enhanced with embedded Process Analytical Technology (ePAT).



### HMI

The Human Machine Interface (HMI) enables intuitive operation thanks to intelligent software support.

### FE35

The specialist for small batches – with the fastest product changeover time in its class.

### FE55

The all-rounder that tackles your tableting challenges: Productive, efficient, and flexible.

### FE75

The double rotary tablet press delivers peak performance for high-volume production of demanding products.

## FE35

The FE35 combines the advantages of the FE Series with the shortest product changeover time in its class. Users benefit from maximum machine availability and investment security.



Tablet		
Number of layers		1-layer
max. tablet output		252,000 units/h
max. tablet output Pmax®		367,200 units/h
max. tablet diameter		25 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	325	325	325	325	325	325	325
<b>Compression force pre- and main compression 1*</b>	kN	0–80	0–80	0–80	0–80	0–80	0–80	0–80
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120	5–120	5–120	5–120	5–120
<b>Dimensions L x W x H</b>	mm	1,026 x 1,042 x 2,043 without integrated switch cabinet 1,336 x 1,042 x 2,043 with integrated switch cabinet						
<b>Weight</b>		Tablet press approx. 2,800–3,000 kg, operating terminal 100 kg, switch cabinet 350 kg						
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, 50/60Hz, power consumption 14 kW						

### Turret specifications

<b>Segments (S) / Dies (D)</b>		S	S	S	S	D	D	D
<b>Number of punch stations</b>		51	33	27	24	35	28	23
<b>Punch type</b>		FS12®	FS19® EU19 TSM19	EU1" TSM1"	EU1"-441	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D
<b>Tablet output units/h</b>		min.	45,900	29,700	24,300	21,600	31,500	25,200
		max.	367,200	237,600	145,800	129,600	252,000	201,600
<b>Max. compression force 1***</b>	kN	34	80	80	80	80	80	80
<b>Max. compression force 2***</b>	kN	34	80	80	80	80	80	80
<b>Max. tablet diameter</b>	mm	11	18	25	25	13	18	25
<b>Max. filling depth****</b>	mm	20	20	20	20	18	18	18
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15
		max.	min <sup>-1</sup>	120	120	90	90	120
<b>Segment / Die height</b>	mm	25	25	25	25	22.225	22.225	23.8
<b>Die diameter</b>	mm					24	30.16	38.1
<b>Punch shaft diameter</b>	mm	12	19	25.35	25.35	19	19	25.35
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6
<b>Upper/lower punch</b>			(133.35)	(133.35)		(133.35)	(133.35)	(133.35)
<b>Upper punch insertion depth</b>	mm	1–5	1–5	1–5	1–5	1–5	1–5	1–5

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request



Brochure: FE Series

# FE55

The FE55 sets benchmarks in three key factors for tableting efficiency: productivity, flexibility, and availability. As standard, it can produce more than 90 percent of all products without additional investments or complex changeovers – including bi-layer tablets.



Tablet	
Number of layers	1-, 2-layer
max. tablet output	417,600 units/h
max. tablet output Pmax®	626,400 units/h
max. tablet diameter	27 mm

## Technical specifications

<b>Pitch circle diameter</b>	mm	550	550	550	550	550	550	550
<b>Compression force pre- and main compression 1*</b>	kN	0–100	0–100	0–100	0–100	0–100	0–100	0–100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120	5–120	5–120	5–120	5–120
<b>Dimensions L x W x H</b>	mm	1,306 x 1,306 x 2,048 without integrated switch cabinet 1,306 x 1,626 x 2,048 with integrated switch cabinet						
<b>Weight</b>		Tablet press approx. 3,700–3,900 kg, operating terminal 100 kg, switch cabinet 350 kg						
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, 50/60 Hz, power consumption 16 kW						

## Turret specifications

<b>Segments (S) / Dies (D)</b>		S	S	S	S	D	D	D		
<b>Number of punch stations</b>		87	60	45	39	58	47	39		
<b>Punch type</b>		FS12®	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D		
<b>Tablet output units/h</b>		min.	78,300	54,000	40,500	–	52,200	42,300	35,100	
		max.	626,400	432,000	243,000	243,000	417,600	338,400	210,600	
<b>Max. compression force 1***</b>	kN	34	100	100	80	100	100	100		
<b>Max. compression force 2***</b>	kN	34	100	100	80	100	100	100		
<b>Max. compression force 3***</b>	kN	34	100	100	–	100	100	100		
<b>Max. tablet diameter</b>	mm	11	18	25	27	13	18	25		
<b>Max. filling depth</b>		1st layer****	mm	20	20	20	20	13	18	18
		2nd layer	mm	8	8	8	–	8	8	8
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15	15	
		max.	min <sup>-1</sup>	120	120	90	90	120	120	90
<b>Segment / Die height</b>	mm	25	25	25	25	22.225	22.225	23.8		
<b>Die diameter</b>	mm					22	30.16	38.1		
<b>Punch shaft diameter</b>	mm	12	19	25.35	28	19	19	25.35		
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6 (		
<i>Upper/lower punch</i>			(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	133.35)		
<b>Upper punch insertion depth</b>	mm	1–5 (8****)	1–5 (8****)	1–5 (8****)	1–5	1–5 (8****)	1–5 (8****)	1–5 (8****)		

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# FE75

The FE75 delivers peak performance for high-volume production of up to 1,656,000 tablets per hour.



Tablet	
Number of layers	1-, 2-layer
max. tablet output	1,080,000 units/h
max. tablet output Pmax®	1,656,000 units/h
max. tablet diameter	27 mm

## Technical specifications

<b>Pitch circle diameter</b>	mm	710	710	710	710	710	710	710
<b>Compression force pre- and main compression 1*</b>	kN	0–100	0–100	0–100	0–100	0–100	0–100	0–100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120	5–120	5–120	5–120	5–120
<b>Dimensions L x W x H</b>	mm	1,463 x 1,463 x 2,046 without integrated switch cabinet 1,463 x 1,778 x 2,046 with integrated switch cabinet						
<b>Weight</b>		Tablet press approx. 5,300–5,500 kg, operating terminal 100 kg, switch cabinet 350 kg						
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, 50/60Hz, power consumption 16 kW						

## Turret specifications

<b>Segments (S) / Dies (D)</b>		S	S	S	S	D	D	D		
<b>Number of punch stations</b>		115	75	55	55	75	63	51		
<b>Punch type</b>		FS12®	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D		
<b>Tablet output units/h</b>		min.	207,000	135,000	99,000	–	135,000	113,400	91,800	
		max.	1,656,000	1,080,000	594,000	594,000	1,080,000	907,200	550,800	
<b>Max. compression force 1***</b>	kN	34	100	100	80	100	100	100		
<b>Max. compression force 2***</b>	kN	34	100	100	80	100	100	100		
<b>Max. compression force 3***</b>	kN	34	100	100	–	100	100	100		
<b>Max. compression force 4***</b>	kN	34	100	100	–	100	100	100		
<b>Max. tablet diameter</b>	mm	11	18	25	27	13	18	25		
<b>Max. filling depth</b>		1st layer****	mm	20	20	20	20	18	18	
		2nd layer	mm	8	8	8	–	8	8	
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15	15	
		max.	min <sup>-1</sup>	120	120	90	90	120	120	90
<b>Segment / Die height</b>	mm	25	25	25	25	22.225	30.16	23.8		
<b>Die diameter</b>	mm					22	30.16	38.1		
<b>Punch shaft diameter</b>	mm	12	19	25.35	28	19	19	25.35		
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.35	133.6	133.6		
<i>Upper/lower punch</i>			(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)		
<b>Upper punch insertion depth</b>	mm	1–5 (8****)	1–5 (8****)	1–5 (8****)	1–5 (8****)	1–5 (8****)	1–5 (8****)	1–5 (8****)		

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# i Series

## Versatile options, no compromises

The i Series from Fette Compacting carries on the tradition of proven, top-quality tablet presses. It offers solutions for higher performance and easy handling – even in challenging productions. Thanks to cross-generation system compatibility and intuitive operation, the i Series simplifies manufacturing and boosts productivity – from start-ups to large pharmaceutical companies.



### HMI

The Human Machine Interface (HMI) enables intuitive operation thanks to intelligent software support.

### F<sub>10i</sub>

The powerful, quickly adaptable single rotary tablet press for smaller batches.

### F<sub>20i</sub>

The all-rounder for a broad spectrum of products and batch sizes.

### F<sub>30i</sub>

The high-performance double rotary tablet press for maximum efficient production of large batches.

## F10i

The F10i is a powerful single rotary tablet press for small batches. Features such as a mechanical-manual turret clamping system ensure flexibility and versatility. This makes the F10i ideally suited for the demands of dynamic production.



Tablet	
Number of layers	1-layer
max. tablet output	230,400 units/h
max. tablet output Pmax®	324,000 units/h
max. tablet diameter	27 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	280	280	280	280	280	280	280	280	280
<b>Compression force pre- and main compression 1*</b>	kN	0–80	0–80	0–80	0–80	0–80	0–80	0–80	0–80	0–80
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120	5–120	5–120	5–120	5–120	5–120	5–120
<b>Dimensions L x W x H</b>	mm	925 x 1,112 x 1,875 ***								
<b>Weight</b>		Tablet press approx. 2,100 kg ****, operating terminal 90 kg								
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, 50/60 Hz, power consumption 9 kW								

### Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S	S
<b>Number of punch stations</b>		32	30	24	20	45	30	24	21	21
<b>Punch type</b>		EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12®	FS19® /EU19 TSM19	EU1" TSM1"	EU1"-441	EU28
<b>Tablet output units/h</b>	min.	48,000	45,000	36,000	30,000	67,500	45,000	36,000	31,500	–
	max.	230,400	216,000	172,800	120,000	324,000	216,000	144,000	126,000	126,000
<b>Max. compression force 1*****</b>	kN	80	80	80	80	34	80	80	80	80
<b>Max. compression force 2*****</b>	kN	80	80	80	80	34	80	80	80	80
<b>Max. tablet diameter</b>	mm	11	13	18	25	11	18	25	25	27
<b>Max. filling depth*****</b>	mm	18	18	18	20	20	20	20	20	20
<b>Recommended turret rotation speed</b>	min.	min <sup>-1</sup>	25	25	25	25	25	25	25	15
	max.	min <sup>-1</sup>	120	120	120	100	120	120	100	100
<b>Die diameter</b>	mm	22	24	30.16	38.1	–	–	–	–	–
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25	25
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	25.35	28
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6
<b>Upper/lower punch</b>		(133.35)	(133.35)	(133.35)	(133.35)		(133.35)	(133.35)		
<b>Upper punch insertion depth</b>	mm	1–5	1–5	1–5	1–5	1–5	1–5	1–5	1–5	1–5

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Dimensions may vary depending on option. \*\*\*\* Weight may vary depending on option. \*\*\*\*\* Limited by punch properties. \*\*\*\*\* Special filling depth available on request



Brochure: i Series

# F20i

With an output of up to 475,200 tablets per hour and the optional Pmax® turret technology, the F20i is the all-rounder of the i Series. It handles a broad spectrum of products and batch sizes, delivering maximum flexibility in manufacturing.



Tablet	
Number of layers	1-layer
max. tablet output	338,400 units/h
max. tablet output Pmax®	475,200 units/h
max. tablet diameter	27 mm

## Technical specifications

<b>Pitch circle diameter</b>	mm	410	410	410	410	410	410	410	410	410
<b>Compression force pre- and main compression 1*</b>	kN	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	1,220 x 1,220 x 2,022								
<b>Weight</b>		Tablet press approx. 3,300-3,500 kg, operating terminal 90 kg, switch cabinet 270 kg								
<b>Electrical supply parameters</b>		Operating voltage 400-480 V, 50/60 Hz, power consumption 13.2 kW								

## Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S	S	
<b>Number of punch stations</b>		47	43	36	30	66	45	36	33	33	
<b>Punch type</b>		EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12® EU19 TSM19	FS19® EU19 TSM19	EU1" TSM1"	EU1"-441 TSM1"	EU28	
<b>Tablet output units/h</b>		min.	42,300	38,700	32,400	27,000	59,400	40,500	32,400	29,700	-
		max.	338,400	309,600	259,200	180,000	475,200	324,000	216,000	198,000	198,000
<b>Max. compression force 1***</b>	kN	100	100	100	100	34	100	100	100	100	80
<b>Max. compression force 2***</b>	kN	100	100	100	100	34	100	100	100	100	80
<b>Max. tablet diameter</b>	mm	11	13	18	25	11	18	25	25	27	
<b>Max. filling depth**</b>	mm	18	18	18	20	20	20	20	20	20	
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15	15	15	15
		max.	min <sup>-1</sup>	120	120	120	100	120	120	100	100
<b>Die diameter</b>	mm	22	24	30.16	38.1	-	-	-	-	-	
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25	25	
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	25.35	28	
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	
<b>Upper/lower punch</b>			(133.35)	(133.35)	(133.35)		(133.35)	(133.35)	(133.35)		
<b>Upper punch insertion depth</b>	mm	1-5	1-5	1-5	1-5	1-5	1-5	1-5	1-5	1-5	

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties

# F30i

The F30i is a high-performance double rotary tablet press with ultra-fast first-layer sampling. Operators can quickly adapt the tablet press to varying requirements. Designed for large-batch production, the F30i delivers a maximum output of up to 1.6 million tablets per hour.



Tablet	
Number of layers	1-, 2-layer
max. tablet output	1,137,600 units/h
max. tablet output Pmax®	1,584,000 units/h
max. tablet diameter	27 mm

## Technical specifications

<b>Pitch circle diameter</b>	mm	680	680	680	680	680	680	680	680	680
<b>Compression force pre- and main compression 1*</b>	kN	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	1,394 x 1,394 x 2,030								
<b>Weight</b>		Tablet press approx. 4,600 kg, operating terminal 90 kg, switch cabinet 350 kg								
<b>Electrical supply parameters</b>		Operating voltage 400-480 V, frequency 50/60 Hz, power consumption 17 kW								

## Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S	
<b>Number of punch stations</b>		79	73	61	49	110	75	55	55	
<b>Punch type</b>		FS19® EU19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12® EU19 TSM19	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28	
<b>Tablet output units/h</b>		min.	142,200	131,400	109,800	88,200	396,000	135,000	99,000	
		max.	1,137,600	1,051,200	878,400	470,400	1,584,000	1,080,000	528,000	528,000
<b>Max. compression force 1***</b>	kN	100	100	100	100	34	100	100	100	80
<b>Max. compression force 2***</b>	kN	100	100	100	100	34	100	100	100	80
<b>Max. compression force 3***</b>	kN	100	100	100	100	34	100	100	100	-
<b>Max. compression force 4***</b>	kN	100	100	100	100	34	100	100	100	-
<b>Max. tablet diameter</b>	mm	11	13	18	20	11	18	25	27	
<b>Max. filling depth****</b>	mm	18	18	18	18	20	20	20	20	
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15	15	15
		max.	min <sup>-1</sup>	120	120	120	80	120	120	100
<b>Die diameter</b>	mm	22	24	30.16	38.1	-	-	-	-	
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25	
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	25.35	28
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	
<b>Upper/lower punch</b>			(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	
<b>Upper punch insertion depth</b>	mm	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# p Series

## Simply efficient production

Whether small batches or high-speed production: The p Series has the right solution for every requirement. You benefit from a uniform technical standard that simplifies your tablet production.



### HMI

Human Machine Interface (HMI): Full control through intuitive touchscreen operation on all models.

### F<sub>10p</sub>

The specialist for precision and fast changeovers with small batches.

### F<sub>20p</sub>

The all-rounder for maximum flexibility across variable batch sizes.

### F<sub>30p</sub>

The powerhouse for high-speed large-scale production.

# F10p

The F10p brings large-scale production efficiency down to small-batch scale. As a single rotary tablet press, it is the ideal choice for development, clinical samples, or niche products. Thanks to identical process technology, scale-up within the p Series is reliable and easy.



Tablet	Number of layers	1-layer
	max. tablet output	230,400 units/h
	max. tablet diameter	27 mm

#### Technical specifications

<b>Pitch circle diameter</b>	mm	280	280	280	280	280	280	280	280
<b>Compression force pre- and main compression 1*</b>	kN	0-80	0-80	0-80	0-80	0-80	0-80	0-80	0-80
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	900 x 1,205 x 1,970							
<b>Weight</b>		Tablet press approx. 2,000 kg, operating terminal 100 kg							
<b>Electrical supply parameters</b>		Operating voltage 380-480 V, 50/60 Hz, power consumption 9 kW***							

#### Turret specifications

<b>Dies (D) / Segments (S)</b>		D	D	D	D	S	S	S	S
<b>Number of punch stations</b>		32	30	24	20	30	24	21	21
<b>Punch type</b>		EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS19® EU19 TSM19	EU1" TSM1"	EU1"-441	EU28
<b>Tablet output units/h</b>	min.	48,000	45,000	36,000	30,000	45,000	36,000	31,500	
	max.	230,400	216,000	172,800	120,000	216,000	144,000	126,000	
<b>Max. compression force 1****</b>	kN	80	80	80	80	80	80	80	
<b>Max. compression force 2****</b>	kN	80	80	80	80	80	80	80	
<b>Max. tablet diameter</b>	mm	11	13	18	25	18	25	25	
<b>Max. filling depth*****</b>	mm	18	18	18	20	20	20	20	
<b>Recommended turret rotation speed</b>	min	min <sup>-1</sup> 25	25	25	25	25	25	25	
	max.	min <sup>-1</sup> 120	120	120	100	120	100	100	
<b>Die diameter</b>	mm	22	24	30.16	38.1	-	-	-	
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	19	25.35	25.35	
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	
<b>Upper/lower punch</b>		(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	
<b>Upper punch insertion depth</b>	mm	1-5	1-5	1-5	1-5	1-5	1-5	1-5	


Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* 14.7 kW maximum power with integrated process equipment. \*\*\*\* Limited by punch properties. \*\*\*\*\* Special filling depth available on request



Brochure: p Series

# F20p

The F20p is the most versatile model of the p Series. As a single rotary tablet press, it covers the broadest spectrum of products and batch sizes. It combines high output with technical sophistication for particularly flexible production – even with demanding formulations.



Tablet	Number of layers	1-layer
	max. tablet output	338,400 units/h
	max. tablet diameter	27 mm

## Technical specifications

<b>Pitch circle diameter</b>	mm	410	410	410	410	410	410	410	410
<b>Compression force pre- and main compression 1*</b>	kN	0–100	0–100	0–100	0–100	0–100	0–100	0–100	0–100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120	5–120	5–120	5–120	5–120	5–120
<b>Dimensions L x W x H</b>	mm	1,200 x 1,545 x 2,080							
<b>Weight</b>		Tablet press approx. 3,300 kg, operating terminal 100 kg							
<b>Electrical supply parameters</b>		Operating voltage 380–480 V, 50/60 Hz, power consumption 12.5 kW ***							

## Turret specifications

<b>Dies (D) / Segments (S)</b>		D	D	D	D	S	S	S	S
<b>Number of punch stations</b>		47	43	36	30	45	36	33	33
<b>Punch type</b>		EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS19® EU19 TSM19	EU1" TSM1"	EU1"-441	EU28
<b>Tablet output units/h</b>	min.	42,300	38,700	32,400	27,000	40,500	32,400	29,700	
	max.	338,400	309,600	259,200	180,000	324,000	216,000	198,000	158,400
<b>Max. compression force 1****</b>	kN	100	100	100	100	100	100	100	80
<b>Max. compression force 2****</b>	kN	100	100	100	100	100	100	100	80
<b>Max. tablet diameter</b>	mm	11	13	18	25	18	25	25	27
<b>Max. filling depth*****</b>	mm	18	18	18	20	20	20	20	20
<b>Recommended turret rotation speed</b>	min.	min <sup>-1</sup> 15	15	15	15	15	15	15	20
	max.	min <sup>-1</sup> 120	120	120	100	120	100	100	80
<b>Die diameter</b>	mm	22	24	30.16	38.1	–	–	–	–
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	19	25.35	25.35	28
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6
<b>Upper/lower punch</b>			(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	
<b>Upper punch insertion depth</b>	mm	1–5	1–5	1–5	1–5	1–5	1–5	1–5	1–5

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* 17 kW maximum power with integrated process equipment. \*\*\*\* Limited by punch properties. \*\*\*\*\* Special filling depth available on request

# F30p

The F30p is the double rotary tablet press of the p Series, purpose-built for large-batch production. It delivers maximum volumes at top speed and offers specialized capabilities for double-layer tablet manufacturing.



Tablet	Number of layers	1-, 2-layer
	max. tablet output	1,137,600 units/h
	max. tablet diameter	27 mm

## Technical specifications

<b>Pitch circle diameter</b>	mm	680	680	680	680	680	680	680	680
<b>Compression force pre- and main compression 1*</b>	kN	0–100	0–100	0–100	0–100	0–100	0–100	0–100	0–100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120	5–120	5–120	5–120	5–120	5–120
<b>Dimensions L x W x H</b>	mm	1,394 x 1,394 x 2,030							
<b>Weight</b>		Tablet press approx. 4,600 kg, operating terminal 90 kg, switch cabinet 350 kg							
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, frequency 50/60 Hz, power consumption 17 kW							

## Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S
<b>Number of punch stations</b>		79	73	61	49	110	75	55	55
<b>Punch type</b>		FS19® EU19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12®	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28
<b>Tablet output units/h</b>	min.	142,200	131,400	109,800	88,200	396,000	135,000	99,000	
	max.	1,137,600	1,051,200	878,400	470,400	1,584,000	1,080,000	528,000	528,000
<b>Max. compression force 1***</b>	kN	100	100	100	100	34	100	100	80
<b>Max. compression force 2***</b>	kN	100	100	100	100	34	100	100	80
<b>Max. compression force 3***</b>	kN	100	100	100	100	34	100	100	–
<b>Max. compression force 4***</b>	kN	100	100	100	100	34	100	100	–
<b>Max. tablet diameter</b>	mm	11	13	18	20	11	18	25	27
<b>Max. filling depth****</b>	mm	18	18	18	18	20	20	20	20
<b>Recommended turret rotation speed</b>	min.	min <sup>-1</sup> 15	15	15	15	15	15	15	15
	max.	min <sup>-1</sup> 120	120	120	80	120	120	100	80
<b>Die diameter</b>	mm	22	24	30.16	38.1	–	–	–	–
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	28
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6
<b>Upper/lower punch</b>			(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	
<b>Upper punch insertion depth</b>	mm	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# Containment Guard

## Processing active ingredients safely

When processing active and highly active ingredients in tablet production, reliable containment is essential to protect operators, the environment, and products. The Containment Guard (CG) from Fette Compacting provides a structured approach to selecting the right Containment solution for your process. Six levels systematically

classify Containment Performance – from low-dust standard production to fully encapsulated high-containment systems. Performance is tested, certified, and documented using a standardized multi-stage exposure measurement procedure. This provides you with a solid foundation even before the final risk assessment.



### CG Level 1–2: Basic Containment

For moderate exposure protection requirements. Closed process chambers and controlled product flow limit particle release. Cleaning and maintenance are performed primarily manually or with simple assisted procedures.

### CG Level 3–4: Enhanced Containment


For elevated protection requirements when handling active ingredients. Improved sealing, secure interfaces, and integrated transfer systems significantly reduce dust emissions. Cleaning is semi-automated or uses prepared systems.

### CG Level 5–6: High Containment

For highly active ingredients with very low exposure limits. Fully encapsulated process chambers, isolator technology, and the patented air management system minimize emissions. Automated wash-in-place systems ensure safe cleaning without manual intervention.

# FE55

## with Containment Package



### Tablet

Number of layers 1-, 2-layer  
 max. tablet output 417,600 units/h  
 max. tablet output Pmax® 626,400 units/h  
 max. tablet diameter 27 mm

### Technical specifications

Pitch circle diameter	mm	550	550	550	550	550	550
Compression force pre- and main compression 1*	kN	0–100	0–100	0–100	0–100	0–100	0–100
Turret rotation speed 2**	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120
Dimensions L x W x H	mm	1,306 x 1,306 x 2,048 without integrated switch cabinet 1,306 x 1,626 x 2,048 with integrated switch cabinet					
Weight		Tablet press approx. 3,700–3,900 kg, operating terminal 100 kg, switch cabinet 350 kg					
Electrical supply parameters		Operating voltage 400–480 V, 50/60 Hz, power consumption 16 kW					

### Turret specifications

Segments (S) / Dies (D)		S	S	S	S	D	D	D
Number of punch stations		87	60	45	39	58	47	39
Punch type		FS12®	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D
Tablet output units/h	min.	78,300	54,000	40,500	–	52,200	42,300	35,100
	max.	626,400	432,000	243,000	243,000	417,600	338,400	210,600
Max. compression force 1***	kN	34	100	100	80	100	100	100
Max. compression force 2***	kN	34	100	100	80	100	100	100
Max. compression force 3***	kN	34	100	100	–	100	100	100
Max. tablet diameter	mm	11	18	25	27	13	18	25
Max. filling depth	1st layer****	mm	20	20	20	20	13	18
	2nd layer	mm	8	8	8	–	8	8
Recommended turret rotation speed	min.	min <sup>-1</sup>	15	15	15	15	15	15
	max.	min <sup>-1</sup>	120	120	90	90	120	90
Segment / Die height	mm	25	25	25	25	22.225	22.225	23.8
Die diameter	mm					22	30.16	38.1
Punch shaft diameter	mm	12	19	25.35	28	19	19	25.35
Punch length	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6 (
Upper/lower punch			(133.35)	(133.35)		(133.35)	(133.35)	133.35)
Upper punch insertion depth	mm	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)	1–5	1–5 (8*****)	1–5 (8*****)	1–5 (8*****)

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation




Brochure: Containment Guard

# FE75

## with Containment Package



	<b>Tablet</b>	
	Number of layers	1-, 2-layer
	max. tablet output	1,080,000 units/h
	max. tablet output Pmax®	1,656,000 units/h
	max. tablet diameter	27 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	710	710	710	710	710	710
<b>Compression force pre- and main compression 1*</b>	kN	0-100	0-100	0-100	0-100	0-100	0-100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	1,463 x 1,463 x 2,046 without integrated switch cabinet 1,463 x 1,778 x 2,046 with integrated switch cabinet					
<b>Weight</b>		Tablet press approx. 5,300-5,500 kg, operating terminal 100 kg, switch cabinet 350 kg					
<b>Electrical supply parameters</b>		Operating voltage 400-480 V, 50/60Hz, power consumption 16 kW					

### Turret specifications


<b>Segments (S) / Dies (D)</b>		S	S	S	S	D	D	D
<b>Number of punch stations</b>		115	75	55	55	75	63	51
<b>Punch type</b>		FS12®	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28	FS19® EU19 TSM19 BB	FS19® EU19 B	EU1" EU1"-441 TSM1" D
<b>Tablet output units/h</b>	min.	207,000	135,000	99,000	-	135,000	113,400	91,800
	max.	1,656,000	1,080,000	594,000	594,000	1,080,000	907,200	550,800
<b>Max. compression force 1***</b>	kN	34	100	100	80	100	100	100
<b>Max. compression force 2***</b>	kN	34	100	100	80	100	100	100
<b>Max. compression force 3***</b>	kN	34	100	100	-	100	100	100
<b>Max. compression force 4***</b>	kN	34	100	100	-	100	100	100
<b>Max. tablet diameter</b>	mm	11	18	25	27	13	18	25
<b>Max. filling depth</b>	1st layer****	mm	20	20	20	18	18	18
	2nd layer	mm	8	8	8	8	8	8
<b>Recommended turret rotation speed</b>	min.	min <sup>-1</sup>	15	15	15	15	15	15
	max.	min <sup>-1</sup>	120	120	90	90	120	90
<b>Segment / Die height</b>	mm	25	25	25	25	22.225	30.16	23.8
<b>Die diameter</b>	mm					22	30.16	38.1
<b>Punch shaft diameter</b>	mm	12	19	25.35	28	19	19	25.35
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.35	133.6	133.6
	Upper/lower punch		(133.35)	(133.35)		(133.35)	(133.35)	(133.35)
<b>Upper punch insertion depth</b>	mm	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)			

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# F10i

## Containment Guard



	<b>Tablet</b>	
	Number of layers	1-layer
	max. tablet output	230,400 units/h
	max. tablet output Pmax®	324,000 units/h
	max. tablet diameter	27 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	280	280	280	280	280	280	280	280	280
<b>Compression force pre- and main compression 1*</b>	kN	0-80	0-80	0-80	0-80	0-80	0-80	0-80	0-80	0-80
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	925 x 1,112 x 1,875 ***								
<b>Weight</b>		Tablet press approx. 2,100 kg****, operating terminal 90 kg								
<b>Electrical supply parameters</b>		Operating voltage 400-480 V, 50/60 Hz, power consumption 9 kW								


### Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S	S
<b>Number of punch stations</b>		32	30	24	20	45	30	24	21	21
<b>Punch type</b>		EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12®	FS19® /EU19 TSM19	EU1" TSM1"	EU1"-441	EU28
<b>Tablet output units/h</b>	min.	48,000	45,000	36,000	30,000	67,500	45,000	36,000	31,500	-
	max.	230,400	216,000	172,800	120,000	324,000	216,000	144,000	126,000	126,000
<b>Max. compression force 1*****</b>	kN	80	80	80	80	34	80	80	80	80
<b>Max. compression force 2*****</b>	kN	80	80	80	80	34	80	80	80	80
<b>Max. tablet diameter</b>	mm	11	13	18	25	11	18	25	25	27
<b>Max. filling depth*****</b>	mm	18	18	18	20	20	20	20	20	20
<b>Recommended turret rotation speed</b>	min.	min <sup>-1</sup>	25	25	25	25	25	25	25	15
	max.	min <sup>-1</sup>	120	120	120	100	120	120	100	100
<b>Die diameter</b>	mm	22	24	30.16	38.1	-	-	-	-	-
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25	25
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	25.35	28
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6
	Upper/lower punch		(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)	(133.35)
<b>Upper punch insertion depth</b>	mm	1-5	1-5	1-5	1-5	1-5	1-5	1-5	1-5	1-5

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Dimensions may vary depending on option. \*\*\*\* Weight may vary depending on option. \*\*\*\*\* Limited by punch properties. \*\*\*\*\* Special filling depth available on request

# F20i

## Containment Guard



Tablet	
Number of layers	1-layer
max. tablet output	338,400 units/h
max. tablet output Pmax®	475,200 units/h
max. tablet diameter	27 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	410	410	410	410	410	410	410	410	410
<b>Compression force pre- and main compression 1*</b>	kN	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	1,220 x 1,220 x 2,022								
<b>Weight</b>		Tablet press approx. 3,300-3,500 kg, operating terminal 90 kg, switch cabinet 270 kg								
<b>Electrical supply parameters</b>		Operating voltage 400-480 V, 50/60 Hz, power consumption 13.2 kW								


### Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S	S	
<b>Number of punch stations</b>		47	43	36	30	66	45	36	33	33	
<b>Punch type</b>		EU19 TSM19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12® EU19 TSM19	FS19® EU19 TSM19	EU1" TSM1"	EU1"-441 EU28	EU28	
<b>Tablet output units/h</b>		min.	42,300	38,700	32,400	27,000	59,400	40,500	32,400	29,700	-
		max.	338,400	309,600	259,200	180,000	475,200	324,000	216,000	198,000	198,000
<b>Max. compression force 1***</b>	kN	100	100	100	100	34	100	100	100	100	80
<b>Max. compression force 2***</b>	kN	100	100	100	100	34	100	100	100	100	80
<b>Max. tablet diameter</b>	mm	11	13	18	25	11	18	25	25	27	
<b>Max. filling depth**</b>	mm	18	18	18	20	20	20	20	20	20	
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15	15	15	15
		max.	min <sup>-1</sup>	120	120	120	100	120	120	100	100
<b>Die diameter</b>	mm	22	24	30.16	38.1	-	-	-	-	-	
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25	25	
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	25.35	28	
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	
<b>Upper/lower punch</b>			(133.35)	(133.35)	(133.35)		(133.35)	(133.35)	(133.35)		
<b>Upper punch insertion depth</b>	mm	1-5	1-5	1-5	1-5	1-5	1-5	1-5	1-5	1-5	

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties

# F30i

## Containment Guard



Tablet	
Number of layers	1-, 2-layer
max. tablet output	1,137,600 units/h
max. tablet output Pmax®	1,584,000 units/h
max. tablet diameter	27 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	680	680	680	680	680	680	680	680	680
<b>Compression force pre- and main compression 1*</b>	kN	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100	0-100
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120	5-120
<b>Dimensions L x W x H</b>	mm	1,394 x 1,394 x 2,030								
<b>Weight</b>		Tablet press approx. 4,600 kg, operating terminal 90 kg, switch cabinet 350 kg								
<b>Electrical supply parameters</b>		Operating voltage 400-480 V, frequency 50/60 Hz, power consumption 17 kW								

### Turret specifications

<b>Die (D) / Segments (S)</b>		D	D	D	D	S	S	S	S	
<b>Number of punch stations</b>		79	73	61	49	110	75	55	55	
<b>Punch type</b>		FS19® EU19 BBS	FS19® EU19 TSM19 BB	FS19® EU19 TSM19 B	EU1" EU1"-441 TSM1" D	FS12® EU19 TSM19	FS19® EU19 TSM19	EU1" EU1"-441 TSM1"	EU28	EU28
<b>Tablet output units/h</b>		min.	142,200	131,400	109,800	88,200	396,000	135,000	99,000	
		max.	1,137,600	1,051,200	878,400	470,400	1,584,000	1,080,000	528,000	528,000
<b>Max. compression force 1***</b>	kN	100	100	100	100	34	100	100	100	80
<b>Max. compression force 2***</b>	kN	100	100	100	100	34	100	100	100	80
<b>Max. compression force 3***</b>	kN	100	100	100	100	34	100	100	100	-
<b>Max. compression force 4***</b>	kN	100	100	100	100	34	100	100	100	-
<b>Max. tablet diameter</b>	mm	11	13	18	20	11	18	25	27	
<b>Max. filling depth****</b>	mm	18	18	18	18	20	20	20	20	
<b>Recommended turret rotation speed</b>		min.	min <sup>-1</sup>	15	15	15	15	15	15	15
		max.	min <sup>-1</sup>	120	120	120	80	120	120	100
<b>Die diameter</b>	mm	22	24	30.16	38.1	-	-	-	-	
<b>Die-/segment height</b>	mm	22.225	22.225	22.225	23.8	25	25	25	25	
<b>Punch shaft diameter</b>	mm	19	19	19	25.35	12	19	25.35	25.35	28
<b>Punch length</b>	mm	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6	133.6
<b>Upper/lower punch</b>			(133.35)	(133.35)	(133.35)		(133.35)	(133.35)	(133.35)	
<b>Upper punch insertion depth</b>	mm	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)	1-5 (8****)

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# Wash in Place

## Clean safely, produce faster



Processing active and highly active substances requires residue-free cleaning between batches. At the same time, cleaning must not compromise productivity. Wash in Place (WiP) combines both: A fully automated, software-controlled cleaning process binds dust, cleans the process chamber emission-free, and cuts downtimes to a minimum. Operators and the environment remain reliably protected at all times.


At the heart of the system is the stationary cleaning unit, seamlessly integrated into the machine control. It precisely doses the cleaning agent, automatically manages the entire cleaning process for the tablet press and peripherals, and ensures that all lines are completely drained upon completion.

The WiP process is available in the 1090i WiP, 2090i WiP, and 3090i WiP tablet presses, covering the full performance spectrum – from smaller batches to flexible high-performance production. Dust-tight systems, contamination-free handling, and integrated cleaning ensure end-to-end containment.

The single rotary press models 1090i WiP and 2090i WiP are designed for flexible production with fast product changeovers. The double rotary press 3090i WiP extends Wash in Place to large batches and is engineered for maximum performance and availability. Integrated glove ports and Rapid Transfer Ports (RTP) provide safe access to the interior. An air management system reliably protects the production area from contamination, even during cleaning.

# 1090i WiP

## Quickly cleaned, quickly changed



### Tablet

Number of layers	1-layer
max. tablet output	216,000 units/h
max. tablet diameter	25 mm

### Technical specifications

<b>Pitch circle diameter</b>	mm	280	280	280
<b>Compression force pre- and main compression 1*</b>	kN	0–80	0–80	0–80
<b>Turret rotation speed 2**</b>	min <sup>-1</sup>	5–120	5–120	5–120
<b>Dimensions L × W × H</b>	mm	960 × 960 × 2,034		
<b>Weight</b>		Tablet press 2,100–2,400 kg, operating terminal 100 kg, switch cabinet 350 kg		
<b>Electrical supply parameters</b>		Operating voltage 400–480 V, 50/60Hz, power consumption 7.7 kW		

### Turret specifications

<b>Segments (S)</b>		S	S	S
<b>Number of punch stations</b>		30	24	21
<b>Punch type</b>		FS19® EU19 TSM19	EU1* TSM1*	EU1*-441
<b>Tablet output units/h</b>	min.	45,000	36,000	31,500
	max.	216,000	144,000	126,000
<b>Max. compression force 1***</b>	kN	80	80	80
<b>Max. compression force 2***</b>	kN	80	80	80
<b>Max. tablet diameter</b>	mm	18	25	25
<b>Max. filling depth****</b>	mm	20	20	20
<b>Recommended turret rotation speed min.</b>	min <sup>-1</sup>	25	25	25
max. (laboratory operation)	min <sup>-1</sup>	120	100	100
<b>Segment height</b>	mm	25	25	25
<b>Punch shaft diameter</b>	mm	19	25.35	25.35
<b>Punch length</b>	mm	133.6	133.6	133.6
Upper/lower punch		(133.35)	(133.35)	
<b>Upper punch insertion depth</b>	mm	1–4	1–4	1–4

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request



Brochure: Containment Guard

# 2090i WiP

## Smart cleaning for higher throughput



Tablet	
Number of layers	1-layer
max. tablet output	324,000 units/h
max. tablet diameter	25 mm

### Technical specifications

Pitch circle diameter	mm	410	410	410
Compression force pre- and main compression 1*	kN	0–100	0–100	0–100
Turret rotation speed 2**	min <sup>-1</sup>	5–120	5–120	5–120
Dimensions L × W × H	mm	1,220 × 1,220 × 2,022		
Weight		Tablet press 3,400–3,600 kg, operating terminal 100 kg, switch cabinet 350 kg		
Electrical supply parameters		Operating voltage 400–480 V, 50/60Hz, power consumption 13 kW		

### Turret specifications

Segments (S)		S	S	S
Number of punch stations		45	36	33
Punch type		FS19® EU19 TSM19	EU1* TSM1*	EU1*-441
Tablet output units/h	min.	40,500	32,400	29,700
	max.	324,000	248,400	198,000
Max. compression force 1***	kN	100	100	100
Max. compression force 2***	kN	100	100	100
Max. tablet diameter	mm	18	25	25
Max. filling depth****	mm	20	20	20
Recommended turret rotation speed min.	min <sup>-1</sup>	15	15	15
max.	min <sup>-1</sup>	120	120	100
Segment height	mm	25	25	25
Punch shaft diameter	mm	19	25.35	25.35
Punch length	mm	133.6	133.6	133.6
Upper/lower punch		(133.5)	(133.35)	(133.35)
Upper punch insertion depth	mm	1–4	1–4	1–4

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request

# 3090i WiP

## Maximum across the board



Tablet	
Number of layers	1-, 2-layer
max. tablet output	1,080,000 units/h
max. tablet diameter	25 mm

### Technical specifications

Pitch circle diameter	mm	680	680	680
Compression force pre- and main compression 1*	kN	0–100	0–100	0–100
Turret rotation speed 2**	min <sup>-1</sup>	5–120	5–120	5–120
Dimensions L × W × H	mm	1,390 × 1,390 × 2,024		
Weight		Tablet press 4,700–5,000 kg, operating terminal 100 kg, switch cabinet 350 kg		
Electrical supply parameters		Operating voltage 400–480 V, 50/60Hz, power consumption 18 kW		


### Turret specifications

Segments (S)		S	S	S
Number of punch stations		75	55	55
Punch type		FS19® EU19 TSM19	EU1* TSM1*	EU1*-441
Tablet output units/h	min.	135,000	99,000	99,000
	max.	1,080,000	528,000	528,000
Max. compression force 1***	kN	100	100	100
Max. compression force 2***	kN	100	100	100
Max. compression force 3***	kN	100	100	100
Max. compression force 4***	kN	100	100	100
Max. tablet diameter	mm	18	25	25
Max. filling depth 1st layer****	mm	20	20	20
2nd layer	mm	8	8	8
Recommended turret rotation speed min.	min <sup>-1</sup>	15	15	15
max.	min <sup>-1</sup>	120	80	80
Segment height	mm	25	25	25
Punch shaft diameter	mm	19	25.35	25.35
Punch length	mm	133.6	133.6	133.6
Upper/lower punch		(133.5)	(133.35)	(133.35)
Upper punch insertion depth	mm	1–4 (8****)	1–4 (8****)	1–4 (8****)

Theoretical values or technical limits: These may vary in practice depending on the product and application. Tablet thickness is a product-dependent variable and can vary significantly. \* The minimum pressing force must exceed a certain value for the punch to rotate the pressure roller. The maximum pressing force depends on the installed tableting tools. \*\* The minimum and maximum possible rotational speed depends on the specific application. \*\*\* Limited by punch properties. \*\*\*\* Special filling depth available on request. \*\*\*\*\* 2-layer-operation

# Tableting Tools

## Always in top form



### Segments

In addition to classic dies in proven quality, Fette Compacting also offers segments. The segment technology replaces the conventional turret die table and delivers significantly higher output with shorter changeover times. Fast changeovers, long service life, and full compatibility with existing production parameters.

- + Faster product changeover through tool-free disassembly
- + Up to 50 percent less product loss
- + Full dimensional accuracy compared to conventional dies

### Punches

Whether standard formats or FS Technology® – our punches are engineered for optimal interplay with the turret and tablet press. High-quality materials and tailored coatings ensure maximum service life.

- + Extended service life and increased dwell time with FS® head geometry
- + Broad coating portfolio to combat sticking, wear, and corrosion
- + Available in widely used international and company-specific standards



### EasyCare

A perfectly coordinated all-in-one solution for transporting, cleaning, and storing your tableting tools. EasyCare enables safe and efficient tool handling while keeping manual interventions to a minimum.

- + TRI.EASY Tool Box System for optimal storage and safe transport
- + GMP-compliant cleaning systems
- + Tool Management System

### Turrets

Turrets are the heart of every tablet press. Our Pmax® turrets with FS12® punches boost your output significantly – without additional machine investments or changes to production parameters.

- + Higher output at consistent tablet quality
- + Individual compression force measurement and individual tablet rejection for seamless quality assurance
- + Proven alternative to multi-tip tooling



Brochure: Tableting Tools

# Process Equipment

Seamless processes in production and active ingredient analysis

## Downstream Processes



**Vertical Deduster**  
Optimal production results through reliable dedusting and gentle deburring – also available in combination with a metal detector and in-process control.



**Metal Detector**  
Maximizing uptime through reliable and fully automatic rejection of metallically contaminated products.



**Gratex**  
Gentle dedusting and deburring directly at the tablet press.

## Process Technology



**PKB (Press Chamber Coating System)**  
Applying powdered lubricant to the press punches during production replaces the need to blend corresponding ingredients into the formulation.

## Containment and WiP



**Isolator**  
Maximum safety for operating personnel when processing highly active products and operating in-process control, dedusting, and metal detection equipment.

## In-Process Control



**4-Parameter Tester – Containment Version**  
The flexible and modular 4-Parameter Tester enables fully automatic measurement of four critical parameters: weight, thickness, diameter, and hardness.

## In-Process Control



**Weightmaster 7.2**  
Automatic weighing system designed specifically for the weight control of tablets and dragees.



**Checkmaster**  
Outstanding tablet quality through fully automated measurement of the key tablet parameters: weight, thickness, diameter, and hardness – with full integration into in-process control.



**AT50**  
Universal tablet testing device for fully integrated in-process control of weight, thickness, diameter, and hardness.



**WT50**  
Compact, robust online weight tester – measures sample weight in real time and automatically adjusts tablet press parameters.

# Handling Systems


## Fast changeovers – keep things moving

### Supplies



#### Handling system

Thanks to electronic drives, a single operator can perform turret changeovers ergonomically and safely.



#### Parts Butler

Space-saving, organized storage of removed parts during maintenance and product changeover cycles enables safe and fast work.



#### Turrets

Fette Compacting turrets stand out for their exceptional quality and can be flexibly deployed for any application.



#### Tool Box

The TRI.EASY Tool Box System comprises three different box types for punches, segments, and dies. In them, tools can be safely stored, transported, and cleaned.




#### Service Carts

The precisely fitting mount ensures secure turret fixation, giving operators full access from all sides for all scheduled work.

# EasyCare


## Safe and efficient handling of tableting tools

Tableting tools can only contribute to ideal tablet quality when transport, cleaning, and storage are handled professionally. To prevent damage and ensure handling is as low-risk and efficient as possible, we offer a comprehensive solution in four steps.



#### Transport

With our all-in-one solution, you cover all tool handling steps safely. The patented TRI.EASY Tool Box System provides precision-fit protection for punches, segments, and dies. A key efficiency feature: Tools are stored in special trays that go directly from the box into the cleaning machine – no time-consuming transfer required.




#### Cleaning

For maximum GMP compliance, the systems from Fette Compacting, ARUNA AG, and Borer Chemie AG work seamlessly together. The TRI.EASY trays with their predefined positions fit precisely into the Aruna cleaning plant. There, the tools are cleaned residue-free with matched cleaning and corrosion protection agents from Borer Chemie.



#### Polishing

Machine polishing with a machine from nortec after cleaning delivers maximum effectiveness and minimizes the risk of damage compared to manual polishing. The food-grade process also prevents tablets from sticking to the punches.



#### Storage

After reconditioning, tools can be stored safely and compactly in the stackable TRI.EASY Tool Box System. They are optimally protected from environmental factors such as moisture and heat. The long-term preservative deconex® HT 1191 extends storage life even further.



Brochure: Tableting Tools

# SmartInterface

Real-time data, anytime  
and anywhere

The server-based SmartInterface gives you access to the production parameters and HMI settings of your tablet press – no matter the time or location. With seamless real-time monitoring on your mobile device, you spot deviations instantly and respond without delay. This way, all production stakeholders maintain full control without needing to be physically at the machine.

#### Your Advantages at a Glance:

- + Maximum responsiveness: No more walking to the machine or suiting up in protective clothing. Check the status remotely and act immediately.
- + Full real-time transparency: An intuitive KPI dashboard and detailed process graphics deliver production and machine status with sub-second latency.
- + Efficient data management: Process data and logs can be recorded during live operation and exported in editable formats without interrupting the machine.
- + Top-tier IT security: Access is provided via your internal network in accordance with your own IT policies. Access rights mirror your standard HMI permissions. No validation-relevant changes are needed for upgrades.



SmartInterface

# Global Training Program

Knowledge that works.  
From lab to production.



Behind every excellent tablet are people who know their machines, understand processes, and make the right decisions.

Our global training program supports you and your team across the entire production life cycle: from formulation development through daily machine operation and maintenance to special applications. Practical, structured, and always focused on unlocking the full potential of your machines and processes.

#### Flexible, Global, Certified

Our experienced trainers work with you and your team either directly at your machines on-site or at one of our worldwide Competence Centers in Brazil, China, Germany, India, and the USA. All training courses are delivered worldwide to uniformly high quality standards and are offered in six languages: English, French, German, Mandarin, Portuguese, and Spanish. Participants complete each course with a proficiency assessment and a certificate according to 21 CFR.

#### Tailored to Your Needs

Every production is unique. That's why we develop customized training programs on request, precisely aligned with your processes, machine configurations, and your team's skill level. Whether adapted content, integration of specific production parameters, or entirely new topics – our experts combine deep technical know-how with instructional expertise, delivering knowledge your team can apply immediately.

#### Your Advantages at a Glance

- + Reduced changeover times and downtime through targeted troubleshooting, changeover, and maintenance
- + Higher output and top standards through optimally configured machines
- + Minimized material losses through solid process and machine knowledge
- + Maximum cost-effectiveness by fully leveraging machines and data



Brochure: Global Training Program

# Next Level Nutrition

## Tableting solutions for food and supplements


No two nutrition formulations are alike. Dietary supplements, moisture-sensitive effervescent tablets, and functional confectionery each place entirely different demands on raw materials and processes. Whether developing the first lab-scale formulation, choosing the right production strategy, or running economical series production: We support you with the right tableting solution, our comprehensive service portfolio, and deep process expertise.

### Your Advantages at a Glance:

- + Process reliability from the start: Our R&D Solutions help you master critical material properties early on.
- + The right machine for every application: Whether small batches, series production, or continuous processing – we offer the optimal tablet press for every product and every requirement.
- + Ultimate flexibility: Easy-to-clean machine concepts based on the Quality by Design principle minimize changeover times during frequent product changes.
- + The right strategy: Together, we define the optimal approach – from flexible batch production to highly efficient continuous manufacturing.
- + Lifetime partnership: We support you throughout the entire production life cycle – from formulation development through scale-up to data-driven process optimization.




### Ready for variety Solutions for every application



**Vitamins and Minerals**  
Dose precisely and produce economically in series production.



**Effervescent and Combination Products**  
Process complex and moisture-sensitive formulations reliably.



**Botanicals**  
Turn naturally varying raw material qualities into robust processes.



**Chewable Tablets, Lozenges, and Confectionery**  
Process demanding textures and sticky blends at high output.




**Multi-Layer and Special Shapes**  
Custom geometries and functionalities for a strong brand presence.



Brochure: Next Level Nutrition





#### Europe

##### **Fette Compacting GmbH**

Grabauer Strasse 24  
 21493 Schwarzenbek, Germany  
 Phone +49 4151 12-0  
[tablet@fette-compacting.com](mailto:tablet@fette-compacting.com)

##### **Fette Compacting Ibérica SL**

Avenida Labradores, 1  
 2ª Planta, Oficina 3  
 28760 Tres Cantos, Spain  
 Phone +34 91 8039689  
[fcib@fette-compacting.com](mailto:fcib@fette-compacting.com)

##### **Fette Compacting France**

10 Allée Bienvenue  
 93160 Noisy-le-Grand, France  
 Phone +33 155 812121  
[fcf@fette-compacting.com](mailto:fcf@fette-compacting.com)

##### **Fette Compacting Italia Srl**

Via Bruno Tosarelli, 155  
 40055 Bologna, Italy  
 Phone +39 340 5603976  
[fcit@fette-compacting.com](mailto:fcit@fette-compacting.com)

##### **Fette Compacting Belgium BV**

Schallenhoevedreef 1b  
 2800 Mechelen, Belgium  
 Phone +32 15 684260  
[fcbe@fette-compacting.com](mailto:fcbe@fette-compacting.com)

##### **EuroPharma Machinery Ltd**

Unit 12 Highview  
 Bordon, Hampshire, GU35 0AX  
 United Kingdom  
 Phone +44 1420 473344  
[admin@europharma.co.uk](mailto:admin@europharma.co.uk)

#### Asia

##### **Fette Compacting (China) Co., Ltd.**

No. 9 Shengtong Road,  
 Moling Sub-District,  
 Jiangning Development Zone,  
 211111 Nanjing  
 Jiangsu Province, P.R.C., China  
 Phone +86 25 52121818  
[fcn@fette-compacting.com](mailto:fcn@fette-compacting.com)

##### **Fette Compacting Machinery India Private Limited**

Plot No. S 115, Phase III B  
 Verna Industrial Estate  
 Verna, Goa 403 722, India  
 Phone +91 90498 54568  
[marketing-in@fette-compacting.com](mailto:marketing-in@fette-compacting.com)

##### **Fette Compacting Asia Pacific Pte Ltd.**

2 Bukit Merah Central #18-01  
 Singapore 159835, Singapore  
 Phone +65 6371 1556  
[infoasiapacific@fette-compacting.com](mailto:infoasiapacific@fette-compacting.com)

##### **Fette Compacting Middle East FZE**

Jebel Ali Free Zone, Jafza  
 Lobby 14, Office 308, Dubai  
 United Arab Emirates  
 Phone +971 4 8808226  
[dubai@fette-compacting.com](mailto:dubai@fette-compacting.com)

#### North America

##### **Fette Compacting America, Inc.**

400 Forge Way  
 Rockaway, NJ 07866, USA  
 Phone +1 973 5868722  
[sales@fetteamerica.com](mailto:sales@fetteamerica.com)

#### Latin America

##### **Fette Compacting America Latina Ltda.**

Av. Cambacica, 610 módulo 10  
 Parque Imperador  
 CEP 13097-160  
 Campinas / SP, Brazil  
 Phone +55 19 37969910  
[fcial-info@fette-compacting.com](mailto:fcial-info@fette-compacting.com)

##### **Fette Compacting Mexico, SA de CV**

Adolfo Prieto No. 1638  
 Colonia Del Valle Sur  
 03100 CDMX, Mexico  
 Phone +52 55 40000653  
[tablet@fette-compacting.com](mailto:tablet@fette-compacting.com)

#### Africa

##### **Fette Compacting North-West Africa**

Résidence les Jardins, Bt B, Chéraga  
 16000 Algiers, Algeria  
 Phone +213 770 117 007  
[fcnwa@fette-compacting.com](mailto:fcnwa@fette-compacting.com)

© by Fette Compacting GmbH

This publication may not be reprinted in whole or part without our express permission.

All right reserved. No rights may be derived from any errors in content or from typographical or typesetting errors. Diagrams, features and dimensions represent the current status on the date of issue of this leaflet. We reserve the right to make technical changes. The visual appearance of the products may not necessarily correspond to the actual appearance in all cases or in every detail.

[www.fette-compacting.com](http://www.fette-compacting.com)