

Pawneer[®]

Eye



A Pet is pure emotion

PAWNEER® EYE is a sterile ophthalmic solution for veterinary use, available in a 5 ml bottle. It is an over-the-counter device for veterinary use.

Application: PAWNEER® EYE is recommended as an adjunct in the treatment of keratitis of traumatic or infectious origin, and in lesions such as corneal ulcers of varying depth (epithelial, stromal, or descemetocele), collagenolytic ulcers, or recurrent corneal epithelial erosions. It is also indicated for maintaining the tear film and lubricating the mucosal surfaces in the management of dry eye and keratoconjunctivitis sicca.

The ingredients contained in PAWNEER® EYE accelerate corneal re-epithelialisation by promoting adhesion between the basal layer of the corneal epithelium and the stroma, and by supporting the realignment of damaged stromal collagen. Thanks to this reorganisation of collagen fibres, the epithelial healing process is accelerated, encouraging the migration of corneal keratinocytes towards the lesion site.

Active Ingredients:

- *Collagen*

PAWNEER® EYE contains $\alpha 1$ chains of recombinant type I collagen protein, which act as collagen precursors by stimulating endogenous synthesis, promoting assembly, and ensuring proper folding. The functionality of collagen supports the physiological restoration of ocular connective tissues. These $\alpha 1$ chains are produced using recombinant DNA genetic engineering technology, ensuring a pure and highly safe ingredient.

- *Acetylcysteine*

Acetylcysteine is well known for its strong mucolytic and antioxidant activity and for restoring the physiological tear film. It facilitates corneal cleansing and helps maintain ocular comfort in keratitis and post-operative conditions. Furthermore, it inhibits metalloproteinase enzymes, including collagenase (which are overexpressed, for example, in herpesvirus infections), thereby promoting proper restoration of ocular tissues.

Key Strengths



In cases of corneal ulcer, **PAWNEER® EYE**, when used in combination with antibiotic therapies, promotes rapid stromal healing and improves corneal transparency.



PAWNEER® EYE has multiple functions with complementary targets, making it suitable for a variety of conditions and enabling the rapid restoration of normal ocular physiology.*



In cases of keratoconjunctivitis sicca (KCS), **PAWNEER® EYE** helps to maintain an optimal level of ocular surface hydration.



Composition: Acetylcysteine, triethanolamine (buffer), methylcellulose, benzalkonium chloride, collagen, excipients q.s.

Directions for Use:

After thoroughly washing your hands, apply 1 or 2 drops into the animal's eye three times a day, or as directed by your veterinary surgeon. Avoid direct contact between the bottle tip and the affected tissue.

Warnings:

PAWNEER® EYE is not a veterinary drug product but a veterinary device indicated for the treatment of ocular wounds or lesions in companion animals such as dogs, cats, and small mammals. The product may cause mild eye irritation or a transient burning sensation.

CONTACTS



+39 02 4775 05 38



info@solutiavetpet.com



SOLUTIA VET & PET S.R.L.
Via Saronnino 70/24,
21040 Origgio (VA)-Italia



solutiavetpet.com

